



SEAL
Analytical



AQ700

**DISCRETE ANALYSER FOR
ENVIRONMENTAL TESTING**



INTELLIGENT DESIGN FOR AUTOMATED ENVIRONMENTAL ANALYSIS

Selecting an analyser to automate your nutrient analysis is less about choosing an analyser with the latest robotics and more about how the analyser design contributes to analysis that is equal or superior to the methods it replaces.

Analysing environmental samples is more demanding than analyzing other sample types and requires a higher priority on very low detection levels, precision and reproducible results. These priorities can only be achieved when an analyser:

- ▶ *Exactly mimics the traditional wet chemistry methods it replaces*
- ▶ *Delivers "equivalency" to approved regulatory methods*
- ▶ *Uses approved spectrophotometric measurement technology*
- ▶ *Includes the accepted and optimum cuvette path length of 10 mm*
- ▶ *Quartz detection cell – optically pure detection technology*
- ▶ *Delivers precision, reproducibility, and low detection limits*
- ▶ *Protects the measurement from signal interference*
- ▶ *Protects samples from cross contamination and carry-over*
- ▶ *Brings the chemical reaction to full completion and steady state*

Combine the right design elements with robotics for automation and you have an analyser that not only automates your wet chemistry, it:

- ▶ *Gives you high speed, high quality data at a low cost per test*
- ▶ *Increases daily volume of samples and range of analytes which can be analyzed*
- ▶ *Significantly reduces cost per test, consuming only microliters of reagents and samples*
- ▶ *Minimises waste disposal costs*
- ▶ *Reduces staff contact with hazardous chemicals*
- ▶ *Substantially reduces the laboratory's overall operating costs*



SEAL
Analytical



HOW DOES A DISCRETE ANALYSER WORK?

A Discrete Analyser completely automates your manual wet chemistry methods, mimicking the operation of a laboratory chemist and adding the ability to measure multiple analytes simultaneously.

A Discrete Analyser will:

- ▶ *Automatically and precisely add sample aliquots and reagent to reaction wells*
- ▶ *Mix*
- ▶ *Wait for the reaction to complete*
- ▶ *Measure the analyte*
- ▶ *Record every step, providing an audit trail*

It should also:


- ▶ *Automatically prepare a calibration from a top standard*
- ▶ *Predilute samples*
- ▶ *Autodilute out-of-range samples*
- ▶ *Autospike samples and report recovery*
- ▶ *Perform true sample blanking*
- ▶ *Automatically insert and run Quality Control (QC) checks*
- ▶ *Link easily with LIMS*

Colourimetric methods can be automated with a Discrete Analyser. With no flow, baseline, peak shapes, pump tubes to monitor, hardware changes or shutdown procedures, your laboratory will achieve true “walk-away” analysis. After a run is finished the Discrete Analyser even washes itself out and enters standby mode.

With miniaturised components the Discrete Analyser needs to use only microliter amounts of reagents and samples, significantly reducing your reagent consumption and associated chemical waste.

SEAL Discrete Analysers will reduce time and errors often associated with manual methods, generate lower cost per test and reduce overall laboratory operating costs and increase efficiency.

SEAL Discrete Analysers are compact, bench top analysers that don't require a fume hood, glassware, pressurisation, cylinder gas, or cooling water, making them the most popular and versatile analysers for environmental labs.





Introducing the SEAL

AQ700



Discrete Analyser

METHODS INCLUDE

Alkalinity
Ammonia
Chloride
Colour
Cyanide
Fluoride
Hardness
Iron
Nitrate/Nitrite
Nitrite
Phenol
Phosphate
Silicate
Sulfate
Total Nitrogen
Total Phosphorous

PLUS MANY MORE

Large sample capacity, high throughput
environmental analyser with four robotic arms enabling
total unattended operation for a long walkaway time –
brings nutrient analysis to a new level.



Analysis according to Standard Methods, EPA,
ASTM, ISO, UKAS and other international standards

BENEFITS INCLUDE

Low detection limits and
excellent reproducibility
using 10 or 20 mm
optical quartz cuvette

Sample size flexibility
from 1.2 mL – 10 mL

Sample capacity from
160 – 480 depending on
vial size. Up to 864 tests.

Total volume per test
only 500 – 600 µL

Cadmium reduction module
for nitrate analysis

Premade reagents available

- ▶ **True unattended operation – including ability to run overnight**
- ▶ **Automated standard preparation and dilution of over-range samples**
- ▶ **Integrated QR / barcode reading for samples and reagents for traceability**
- ▶ **Varying size sample trays available to accommodate different workloads**
- ▶ **Segregated chemical waste and wash minimises waste disposal**
- ▶ **LIMS compatible – export in .csv format**

Designed by chemists **for chemists.**



AQ700

**Highest speed and capacity.
Lower detection levels.**

TESTS / CHEMISTRIES

Simultaneous Chemistries	1 - 20
Tests Programmable Per Sample	YES
Test Capacity	864
Cadmium Coil	Integrated
Total Volume Per Test	500 - 600 µL (sample & reagent)

SAMPLES

Sampling Rate	Subject to chemistry
Sample Blanking	YES
Add Samples After Run Commenced	YES
Sample Trays (removable)	4 (40, 106)
Sample Capacity	160 - 480
Sample Consumption	2 - 500 µL
Sample Cup Sizes	1.2 mL, 2 mL, 5 mL, 10 mL
Sampling Arms	3

REAGENTS

Reagent Capacity	24
Reagent Cooling	YES
Reagent Volume	40 mL
Reagent Wells	Disposable
Reagent Monitoring	Automatic
Reagent Level Sensing	YES

OPERATION

Auto Start-up & Shut-down	YES
Auto-dilution	YES
Automated Spike Preparation	YES
Automated Standard Preparation	YES
Segregated Wash Waste	YES

TECHNOLOGY

Wash Stations	4
Integrated sampling probe washer	YES
Cuvette Path Length	10 mm - 20 mm
Optically Pure Cuvette	Optical Quartz Cuvette
Cuvette Cleaning	Automatic
Barcode Reader	YES (sample & reagent)
Simplified Access For Maintenance	YES
Detector	Stationary measurement cell
Filter Wheel	9 filter positions, 350-880 nm
Lamp	Quartz Tungsten - Halogen

SOFTWARE

Data Output	LIMS compatible. Export in .csv
Software Updates	Free
Requirements	Windows version 7 or later

SPECIFICATIONS

Bench-top Analyser	YES
Dimensions (cm)	120 W x 90 D x 90 H
Weight	115 kg
Power Requirements	110 V 60 HZ or 220-240V 50 HZ. Configurable.

FAST, ON-DEMAND ANALYSIS

*Easy, rapid colourimetric testing
with minimal start-up time.*



INTEGRATED OPTICAL QUARTZ CUVETTE

10 mm pathlength or longer for maximum sensitivity and lower detection levels. Quartz is superior to styrene for sample analysis ensuring highest precision.

LOWER DETECTION LEVELS

Critical for environmental applications, lowest possible detection levels are a priority. This is made possible with the right combination of mixing technique, longer path length, optically pure detection, accurate dispensing and completion of chemical reaction.

NO CROSS CONTAMINATION

The only discrete analyser with integrated probe washer. Eliminates cross contamination between reagents and samples. Keeps the probe free of reagents, oil and grease. Ideal for all water, wastewater and soils.



REAGENT WEDGES

With onboard cooling; built-in level sensor to verify reagent volume required for each test.

EFFECTIVE SAMPLE & REAGENT MIXING

Reproducible results thanks to thorough sample and reagent mixing that approximates manual mixing in a flask.



DISPOSABLE REACTION WELLS

Inexpensive, disposable wells that reduce carryover and cost per test. Recyclable.

INTEGRATED CADMIUM COIL

Allows flexibility in nitrate + nitrite testing. Software automatically switches the coil inline. All 4 x common approved nitrate + nitrite chemistry options available. In-situ regeneration.



COMPLETE REACTION

Constant heating and programmable reaction time for a highly controlled reaction. This means the reaction is brought to completion increasing precision and accuracy of test results.

LOW REAGENT CONSUMPTION & WASTE GENERATION

Uses only µL dispenses of reagents and samples to greatly reduce the amount of chemical used and waste generated with each test.




EXTERIOR SEGREGATED WASTE MANAGEMENT SYSTEM

Segregated chemical and wash waste minimises hazardous waste disposal costs. Easy to access and outside of instrument.

► **LIMS READY** Customizable output for easy integration.

USEPA, ASTM, ISO APPROVED METHODS

Also complies with other international regulatory methods.



Software Designed for Environmental Laboratories

AQ700

User-friendly, intuitive, highly flexible software streamlines run set up

Continuous in-house development incorporates user requested features

Controls all analytical procedures from working standard dilution to sample analysis, cuvette washing and system QC

Automated system quality control with built in QCPro™

User can specify Quality Control types, limits and automated corrective actions upon failure

FEATURES

- ▶ Real-time monitoring of reagents
- ▶ Provides an audit trail of all sample analysis
- ▶ Prepares working standards from a stock solution
- ▶ Prepares spiked samples and calculates recoveries
- ▶ Automatic rerun of over-range samples. Diluted over-range samples will be batched with associated QC needed for reportability
- ▶ Intuitive Range Switching - over or under range sample results matched to other calibration curves intuitively by the software
- ▶ Data exportable to LIMS or worksheets
- ▶ Assigns tests in the highest order of efficiency
- ▶ Automatically performs system calibration and general maintenance
- ▶ Easily monitor run status with colour coding to visually indicate reagent, sample and test status
- ▶ Quickly run multiple tests in any order
- ▶ Predicts when analysis will be completed for better task planning
- ▶ Continuously monitors analyser status and temperature of reaction ring

Colourimetric Nutrient Analysers

DISCRETE ANALYSERS



AQ300



AQ400



AQ700

SEGMENTED FLOW ANALYSERS



AA100



AA500



QuAAtro39

50 Years of Experience in Environmental Analysis Built into Every Analyser

50 years' experience in designing, developing and manufacturing automated wet chemistry analysers specifically for very low detection levels in environmental applications has helped SEAL to apply the most useful, easy to use features into the SEAL range of Discrete and Segmented Flow analysers. The SEAL analysers are widely acknowledged as the best for environmental analysis, giving you everything you need to achieve equal or superior results to the manual and approved laboratory methods the SEAL analyser replaces.

SEAL Analysers are monitoring environmental samples in every corner of the globe. They are manufactured in the USA, Germany and the Netherlands. Engineering and chemistry support is provided from SEAL global facilities in USA, Germany, UK, the Netherlands and China along with a worldwide network of specialist distributors.

COMPREHENSIVE SUPPORT


We offer comprehensive applications, technical service and software support.

INCLUDING

- ▶ A choice of preventative maintenance and service contracts to meet your specific requirements
- ▶ In-house and online training
- ▶ Guaranteed availability of genuine consumables and spare parts
- ▶ Adaptation of methods to specific requirements such as matrix, range or detection limit
- ▶ Continuous in-house development of software to incorporate new customer requested features

Robotic Handling Systems

SEAL Robotic MiniLab systems for automating sample pretreatment in the laboratory — improving your sample handling efficiency. Typical applications include BOD, pH, COD, Alkalinity, and conductivity measurements with options such as decapping/capping, sample splitting, and filtration. Call us about your laboratory needs and we will design a robot to suit you.



SEAL MiniLab

Digestion Systems


FOR METALS AND TKN,
TP DIGESTION



BD50



SmartBlock II



DEENA

SEAL
Analytical



www.seal-analytical.com

SEAL Analytical is a global company with offices worldwide - contact us at:

SEAL Analytical, Inc.
6501 West Donges Bay Road
Mequon, WI 53092
United States
Tel: +1 (262) 241 7900
Fax: +1 (262) 241 7970
sales@seal-us.com

SEAL Analytical Ltd.
Porvair Sciences
Clywedog Road South
Wrexham Industrial Estate
Wrexham LL13 9XS
United Kingdom
Tel: +44 (0) 1978 807273
sales.uk@seal-analytical.com

SEAL Analytical Netherlands
ROHASYS BV
Provincienbaan 4
5121 DL Rijen
The Netherlands
Tel: +31 161 240152
Fax: +31 161 240153
salesnl@seal-analytical.com

SEAL Analytical GmbH
Werkstrasse 5
D-22844 Norderstedt
Germany
Tel: +49 (0)40 60 9292 9-00
Fax: +49 (0)40 60 9292 9-02
salesde@seal-analytical.com

SEAL Analytical Shanghai
Room 413, 12th Building,
No. 128 Xiangyin Road,
Shanghai, 200433
China
Tel: +86 21 3362 5002
Fax: +86 21 3362 5002

AQ700 4.0Q EU 1221