



Personal Protective Equipment

INSTRUCTIONS



1. Wash hands before using
2. Grab both cords in one hand and the front of the mask in the other
3. Pull over your head and position as show below with one cord going over your ears
4. Pull barrel lock so it fits snug
5. Once the the mask feels snug, pinch the inner wire on the nose area to get tight fit (see below)

Be Safe!



WASHING INSTRUCTIONS

Hand wash in Cold Water with detergent or soap and hang dry
Always wash your hands before and after using your PPE mask

DO NOT MACHINE WASH




Personal Protective Equipment

INSTRUCTIONS



1. Wash hands before using
2. Grab both cords in one hand and the front of the mask in the other
3. Pull over your head and position as show below with one cord going over your ears
4. Pull barrel lock so it fits snug
5. Once the the mask feels snug, pinch the inner wire on the nose area to get tight fit (see below)

Be Safe!



WASHING INSTRUCTIONS

Hand wash in Cold Water with detergent or soap and hang dry
Always wash your hands before and after using your PPE mask

DO NOT MACHINE WASH