

HLA-FluoGene kits

ARTICLE NO.	PRODUCT	TESTS/PLATE	TESTS/KIT
002 070 010	HLA-FluoGene A TaqMan® based typing of HLA-A alleles	1	10
002 071 010 / 20	HLA-FluoGene B TaqMan® based typing of HLA-B alleles	1 / 2	10 / 20
002 072 010	HLA-FluoGene C TaqMan® based typing of HLA-C alleles	1	10
002 073 010	HLA-FluoGene ABC TaqMan® based typing of HLA class I alleles	1	10
002 074 010	HLA-FluoGene ABDR TaqMan® based typing of HLA-A, -B, -DRB1, 3, 4 & 5 alleles	1	10
002 075 010	HLA-FluoGene DR TaqMan® based typing of HLA-DRB1, 3, 4 & 5 alleles	1	10
002 077 010 / 30	HLA-FluoGene DRDQ TaqMan® based typing of HLA-DRB1, 3, 4, 5 & -DQB1 alleles	1 / 3	10 / 30
002 081 010	HLA-FluoGene DRDQP plus TaqMan® based typing of HLA-DRB1, 3, 4 & 5, -DQA1/-DQB1 and -DPA1/-DPB1 alleles	1	10
002 082 010	HLA-FluoGene^{NX} ABCDRDQ TaqMan® based typing of HLA-A, -B, -C, -DRB1, 3, 4 & 5, -DQA1 and -DQB1 alleles	1	10
002 083 005	HLA-FluoGene^{NX} Match TaqMan® based typing of HLA-A, -B, -C, -DRB1, 3, 4 & 5, -DQA1, -DQB1, -DPA1 and -DPB1 alleles in SABR resolution. Detection of all CWD (2.0) null alleles and all CIWD (3.0) common and intermediate null alleles.	1*	5
002 085 010 / 20	HLA-FluoGene³⁸⁴ Match TaqMan® based typing of HLA-A, -B, -C, -DRB1, 3, 4 & 5, -DQA1, -DQB1, -DPA1 and -DPB1 alleles in SABR resolution. Detection of all CWD (2.0) null alleles and all CIWD (3.0) common and intermediate null alleles.	1 / 2	10 / 20
002 078 048	HLA-FluoGene B27 TaqMan® based typing of <i>HLA-B*27</i> alleles associated with Morbus Bechterew	12	48
002 078 00Q	HLA-FluoGene B27 Q TaqMan® based typing of <i>HLA-B*27</i> alleles associated with Morbus Bechterew Liquid mix compatible with multiple qPCR machines	12	48

* 1 test on 2 plates

Additional FluoGene items

ARTICLE NO.	PRODUCT	TESTS/PLATE	TESTS/KIT
002 07W T96	FluoGene Wipe Test	12	48
002 7WT 384	FluoGene Wipe Test 384	12	48
006 010 005	FluoGene Color Control	1	10
005 07C O010	FluoGene Color Control 384	1	10

HPA-FluoGene kits

ARTICLE NO.	PRODUCT	TESTS/PLATE	TESTS/KIT
003 015 048 / 10	HPA-FluoGene ^{NX} HPA-1 a/b, -2 a/b, -3 a/b, -4 a/b, -5 a/b, -6 a/b -9 a/b, -15 a/b	6 / 1	48 / 10
003 011 048	HPA-FluoGene 1 a/b Screen	12	48

RBC-FluoGene kits

ARTICLE NO.	PRODUCT	TESTS/PLATE	TESTS/KIT
001 086 040 / 10	RBC-FluoGene ABO basic 01, 02, B, A1, A2	4 / 1	40 / 10
001 080 010	RBC-FluoGene ^{NX} ABO plus 01, 02, B, A1, A2, A variant (e. g. Aw, Ael), B variant (e. g. Bw, Bel), 0 variant	1	10
001 081 030 / 10	RBC-FluoGene vERYfy <i>RHD</i> : exons 3, 5, 10, psi; <i>RHCE</i> : C, C ^w , c, E, e; KEL1(K), KEL2(k), JK1(Jk ^a), JK2(Jk ^b), FY1(Fy ^a), FY2(Fy ^b), FY(a- b-), FYX(Fy ^{bweak}), MNS1(M), MNS2(N), MNS3(S), MNS4(s), U+var(P2), U+var(NY), DO1(Do ^a), DO2(Do ^b)	3 / 1	30 / 10
001 091 010	RBC-FluoGene vERYfy eXtend <i>RHCE</i> : C, C ^w , C ^x , c, E, e, VS+, V-, ceAG; <i>RHD</i> : exons 1-7, 9, 10, psi, DNB, D cat VII, DHMi, DAU, weak D Typ 1, 2, 3, 4.0/ 4.3, 4.1, 4.2.1/4.2.2, 5, 11, 14, 15, 17, DEL1, DEL8; KEL1(K), KEL2(k), KEL3(Kp ^a), KEL4(Kp ^b), KEL6(Js ^a), KEL7(Js ^b), JK1(Jk ^a), JK2(Jk ^b), FY1(Fy ^a), FY2(Fy ^b), FY(a- b-), FYX(Fy ^{bweak}), MNS1(M), MNS2(N), MNS3(S), MNS4(s), MNS5(U+), U+var(P2), U+var(P3), U+var(NY), U- (GYPB*05N.01/*05N.02), DO1(Do ^a), DO2(Do ^b), DO:-4(Hy-), DO:-5(Joa-), LU1(Lu ^a), LU2(Lu ^b), YT1(Yt ^a), YT2(Yt ^b), DI1(Di ^a), DI2(Di ^b), Vel+, Vel-, CO1(Co ^a), CO2(Co ^b), KN1(Kn ^a), KN2(Kn ^b) HBB/HBS	1	10
001 082 040 / 10	RBC-FluoGene CDE <i>RHCE</i> : C, C ^w , c, E, e, 733G, VS+,ceAG; <i>RHD</i> : exons 1-7, 9, 10, psi, DNB, DVII, DHMi, DAU0, non-DAU, DAU1, DAU2/6/7, DAU3/7/11	4 / 1	40 / 10
001 092 010	RBC-FluoGene CDE eXtend <i>RHCE</i> : C, C ^w , C ^x , c, E, e, ceAG, ceJAL, CeRN, ceMO, ceAR, ceSM, ceBI, ceEK, ceTI, VS, V, hrS-, hrB-, Crawford+, JAL+, <i>RHD</i> : RHD*01N.01(d) via Hybridbox , exons 1-7, 9, 10, psi, DNB, DIIIa, DIIIa-CEVS(4-7)-D, DVIIa, DVII, DHMi, DAU0, non-DAU, DAU1, DAU2/6/7, DAU3/7/11, weak D type 1, 2, 3, 4.0/ 4.3, 4.1, 4.2.1/4.2.2, 5, 11, 14, 15, 17, 59, DEL1, DEL8	1	10
001 083 040/ 10	RBC-FluoGene D weak/variant weak D type 1, 2, 3, 4.0/ 4.3, 4.1, 4.2.1/4.2.2, 5, 11, 14, 15, 17, DEL1, DEL8	4 / 1	40 / 10
001 085 010	RBC-FluoGene Rare DI1(Di ^a), DI2(Di ^b), DI3(Wr ^a), DI4(Wr ^b), YT1(Yt ^a), YT2(Yt ^b), LU1(Lu ^a), LU2(Lu ^b), KEL3(Kp ^a), KEL4(Kp ^b), KEL6(Js ^a), KEL7(Js ^b), CO1(Co ^a), CO2(Co ^b), KN1(Kn ^a), KN2(Kn ^b), DO1(Do ^a), DO2(Do ^b), DO:-4(Hy-), DO:-5(Joa-), LW5(LW ^a), LW7(LW ^b)	1	10
001 087 048	RBC-FluoGene D-Screen Screening of <i>RHD</i> exons 3/5 and 10 in a single reaction.	12	48
001 089 010	RBC-FluoGene KKD KEL1(K), KEL2(k), KEL3(Kp ^a), KEL4(Kp ^b), KEL6(Js ^a), KEL7(Js ^b), JK1(Jk ^a), JK2(Jk ^b), FY1(Fy ^a), FY2(Fy ^b), FY(a- b-), FYX(Fy ^{bweak})	1	10



V1.3 QMS 07.25



For USA and Canada:
For Research Use Only.
Not for use in diagnostic procedures.

Contact worldwide

inno-train Diagnostik GmbH
Niederhöchstädter Str. 62
61476 Kronberg
Germany

Phone +49 (0)6173- 6079- 30
Fax +49 (0)6173- 6079- 50
info@inno-train.de

North American Division

inno-train Diagnostik
108 Utica Street
Clinton NY 13323
USA

Phone +1 (888) 352 2196
Fax +1 (315) 381 5609
info@innotrain.org

Linkedin



Homepage

