

# RWJBarnabas HEALTH

## Dragon PowerMic Mobile Configuration

- Access the **Apple app store** or **Google Play store** and Download and Install **PowerMic Mobile**
- Open the app and accept the user agreement.
- Scan the QR Code below, with your smart phone camera.
- Select either **iOS** or **Android** seen at the bottom of the page.



### Getting Started: PowerMic Mobile

Mobile App Configuration

PowerMic Mobile turns your smartphone into a secure wireless microphone for use with Windows-based desktop clinical speech recognition solutions. Optimized for use with Dragon Medical One and other Nuance healthcare solutions, PowerMic Mobile gives clinicians the freedom to roam from workstation to workstation, room to room, and location to location to complete clinical documentation using a smartphone as a wireless microphone.

PowerMic Mobile must be configured to access your institution's organization in the hosted Nuance Management Server.

Complete these steps to configure the PowerMic Mobile app:

1. Provide the address to this page to providers at your site.
2. Have the providers access this page on their mobile devices.
3. Have the providers tap the link below that corresponds to their mobile operating system.

The link will open and configure the PowerMic Mobile app.

PowerMic Mobile Configuration Links

iOS

Android



- Log into the Dragon application from Epic first, then log into the PowerMic Mobile app using your **Dragon Application Login ID** (Same as Epic/Network username). Device should sync (approx. 2-3 secs)
- Launch PMM from your phone and login in using your network username.



(Not Connected)



(PMM Connected)



(Ready to start dictating)

