




purvigor

N99
PROTECTIVE MASK



3D
DESIGN

- FFP3 PROTECTION
- 5 LAYERS OF PROTECTION
- ANTI-INFECTION
- DUST PROTECTION PM 2.5
- This product is a NON-MEDICAL device.

10 PCS


Non-woven Filter, Non-woven Fabric

Attention

Keep in room temperature and do not store with any toxin material.
Do not discard into toilet.

1. Open the face mask, hold the ear straps with both hands and let the nose clip be above the mask.
2. Keep the lower part of the mask close to your chin to cover your mouth and nose.
3. Hang on the ear straps to both ears separately.
4. Press the nose clip with both hands. Press from both sides of the nose to the face until they are completely pressed into the shape of the nose bridge.

1. Open the face mask, hold the ear straps with both hands and let the nose clip be above the mask.
2. Keep the lower part of the mask close to your chin to cover your mouth and nose.
3. Hang on the ear straps to both ears separately.
4. Press the nose clip with both hands. Press from both sides of the nose to the face until they are completely pressed into the shape of the nose bridge.




5-layers
Strengthened Filtration



Application: Helps filter out dust, bacteria,
Manufacturer: GuangDong Fei Fan Mstar
Manufacturer: A1 workshop, No.10,
Jinshui Avenue, zone 2, Caobian No.3
Town, Nanhai District, Foshan
City, Guangdong, China
Date: (dd/mm/yy)

Standard (China)
Standard (European Union)






purvigor

N99

PROTECTIVE MASK



3D

DESIGN

• FFP3 PROTECTION

• 5 LAYERS OF PROTECTION

• ANTI-INFECTION

• DUST PROTECTION PM 2.5

• This product is a NON-MEDICAL device.

1 pc







00:20

飞梵
purvigor

产品合格证

产品名称/Product name: N99口罩(非医用)
Protective Mask N99(Non-Medical Use)

产品型号/Model no: 009

产品规格/Product specification: 15.5 × 10.5cm

检验/Test: 合格Qualified

产品性能/Product performance: PFE ≥ 99%

原材料/Raw materials: 40% 无纺布 + 60% 熔喷布
40% Non-woven fabric + 60% Melt blown fabric

生产批号/Batch number: 2020041502

生产日期/Date of manufacture: 20200425

有效期/Term of validity: 2 years



广东飞梵材料科技有限公司

Guangdong Feifan Material Technology Co., Ltd

地址: 广东省佛山市南海区大沥镇曹边第三工业区
Address: Caobian Third Industrial Zone, Dali Town,
Nanhai District, Foshan City, Guangdong Province