

Exhibit A

Loop Mask (Model – MN)



Model-MN Face Mask

- 2 layers of 100% Polypropylene (Fluid Resistant)
- Breathable
- Soft inner surface
- Adjustable Nose Clip

Donning Instructions	Technical Data								
<ol style="list-style-type: none">1. Hang loop of mask around neck with material resting in front of you.2. Grab one strap with each hand and pull upward to snug the straps.3. Flip the mask upwards covering your face.4. Secure nose clip.	<table border="1"><tbody><tr><td data-bbox="800 1606 992 1682">Material</td><td data-bbox="992 1606 1461 1682">2 Layers 100% Polypropylene 1.25 oz/yd².</td></tr><tr><td data-bbox="800 1682 992 1730">Seams</td><td data-bbox="992 1682 1461 1730">Tex-60 Polyester Thread</td></tr><tr><td data-bbox="800 1730 992 1778">Nose Clip</td><td data-bbox="992 1730 1461 1778">Metal Clip/Plastic Double Wire Tie</td></tr><tr><td data-bbox="800 1778 992 1869">Straps</td><td data-bbox="992 1778 1461 1869">Non-woven Spunbond Polypropylene</td></tr></tbody></table>	Material	2 Layers 100% Polypropylene 1.25 oz/yd ² .	Seams	Tex-60 Polyester Thread	Nose Clip	Metal Clip/Plastic Double Wire Tie	Straps	Non-woven Spunbond Polypropylene
Material	2 Layers 100% Polypropylene 1.25 oz/yd ² .								
Seams	Tex-60 Polyester Thread								
Nose Clip	Metal Clip/Plastic Double Wire Tie								
Straps	Non-woven Spunbond Polypropylene								

This is not a respirator – this is intended to be used as a loop mask.