

- 1** Remove your hearing aids at night.



- 2** Open the battery doors all the way every night.

- 3** Every night, use a clean, dry tissue to gently wipe the outer part of your hearing aids. If you have earmolds, use the cleaning tool loop to clean the earmold's openings.



- 4** Keep hearing aids away from extreme heat and cold.



- 5** Do not wear hearing aids while applying makeup or hair products.



- 6** Keep hearing aids and batteries away from children and pets.



- 7** Replace batteries when low battery beep occurs. Your size 13 battery will last approximately 2 weeks, size 312 battery: approximately 7 days, and size 10 battery: approximately 5 days.



- 8** Remove the sticker from the back of the battery. Place it on your calendar every time you replace your battery. Wait a couple of minutes after removing the sticker before inserting the battery into your hearing aid.



- 9** Please schedule a checkup appointment every 3 months. We will inspect your ears for wax, as well as professionally clean and check your hearing aids to ensure that you continue to receive the full benefits from your devices.



- 10** Get in the healthy habit of daily hearing aid use and enjoy life to the fullest.



IT IS STRONGLY RECOMMENDED YOU USE A HEARING AID DRYER EVERY NIGHT

● Blue = Left Hearing Aid 1 Beep = 2 Beeps = 3 Beeps = ● Red = Right Hearing Aid