

- 1** Remove your hearing aids at night.



- 2** Place devices on charger nightly, making sure both lights activate to indicate charging.




- 3** Nightly, use a clean, dry tissue to gently wipe the outer part of hearing aids, especially bottom section contacts.



- 4** Keep hearing aids away from extreme heat and cold.



- 5** Do not wear hearing aids while applying makeup or hair products.



- 6** Keep hearing aids and batteries away from children and pets.



- 7** Use your charger in a dry place, not near shower or water.



- 8** Daily, remove your hearing aids from charger and clean as instructed before putting them on.



- 9** Please schedule a checkup appointment every 4 months. We will inspect your ears for wax, as well as professionally clean and check your hearing aids to ensure that you continue to receive the full benefits from your devices.



- 10** Get in the healthy habit of daily hearing aid use and enjoy life to the fullest.



● Blue = Left Hearing Aid 1 Beep = 2 Beeps = 3 Beeps = ● Red = Right Hearing Aid