


HELP PREVENT HEARING AID LOSS OR DAMAGE

- 1) Visually check for aids whenever transferring patient to or from the bed.



- 2) Visually check for aids when changing the sheets to prevent washing the hearing aids in the laundry.



- 3) Wait two minutes after peeling sticker off battery before putting it in the hearing aid.

SUN	MON	TUE
1	2	3
8	9	10

- 4) Place sticker on calendar to track usage.

- 5) NEVER place hearing aids in tissues, napkins or paper towels.



AN OUNCE OF PREVENTION
IS WORTH A POUND OF CURE!