

PCSE ONLINE – INFORMATION SHEET

	Optinet Version 5	FLEX	FLEX (RDS)
SUBMISSION METHOD			
PCSE Online	✓	✓	✓
<p><i>Web-based solution for submitting GOS forms online</i></p> <ul style="list-style-type: none"> • Users complete GOS forms online • Allows for signature to be captured electronically • Claims validated in real-time resulting in: <ul style="list-style-type: none"> ◦ Fewer rejected claims ◦ Greater accuracy of payments • Ability to track the status of claims • Online statements are more detailed than current paper versions, which will significantly help the reconciliation process. 			
eGOS	x	✓	✓
<p><i>Solution to allow FLEX to submit GOS claims electronically to PCSE</i></p> <ul style="list-style-type: none"> • Removes the need to double key because patient information will be populated directly from the details already in FLEX • Signatures are captured electronically • Any rejections will be notified within 15 minutes of submission allowing time for details to be corrected and re-submitted before the cut-off date. 			

HARDWARE REQUIRED FOR SIGNING PROCESS (applies to eGOS only)

Patient & Performer



REQUIRED

Signotec Sigma (#ST-BE105-2-FT100)

- 4" (10.5cm) high resolution display
- Backlit to enable easier reading of text
- 2m USB cable
- Price £195.00 approx. (€213.00 Inc. delivery)

Order direct, using the options displayed below:

<https://www.signotec-shop.com/webstore/en/signature-pad-sigma.html>

Product No.:	ST-BE105-2-FT100
Shipping time:	● approx. 1 - 3 working days
Backlight:	with backlight ▼
Connection :	FTDI (virtual COM) ▼
Warranty:	2 years (inclusive) ▼