

Digital Accessories

Be kind to your patients and the environment with new COMFEE BYTES!

۲

Comfee Disposable bite block for 3D imaging

These Single use Bite covers are a comfortable and inexpensive replacement for the Sirona^{*} XIOS foam bite blocks, Sirona item 6141449.

Easy to use, and completely biodegradable, new Comfee Bytes are the environmentally friendly and economical alternative to Sirona's* XiosFoam bite blocks. Comfee Bytes fold over the fixed bite block to form a soft, comfortable bite plate for patients to bite on.

Comfee Bytes will work with almost any 3D CT machine that includes a bite block for patient positioning.

Single use. Affordable. Biodegradable.

Order a box or request a sample today!



Item No. 81499 Comfee Bytes Box of 100 pcs

For ordering information and free samples please visit www.flowdental.com Or call 1-800-356-9729

*Sirona is a registered trademark of Dentsply Sirona.






REPLACE



800.356.9729 | Fax: 631.242.1001 | www.flowdental.com





Safety & Digital Accessories

Other Great CBCT Imaging Products

CBCT Aprons

Available in both Lead free and conventional lead, this apron style is perfect for low dose CBCT X-Rays. The X-Rays from a CBCT travel in parallel waves, so this shorter apron length provides all the protection you need.

The poncho design completely circles your patient's neck and shoulders to provide protection during the full rotation of the X-Ray Tube. Available with or without attached Thyroid collar, in lead or lead free (NL)

75090 and 75090NL for Apron with collar 75091 and 75091NL for Apron without collar

Economically priced and available in 14 colors and 4 attractive patterns

Bite Block Baggies

Avoid costly and time consuming sterilization of your CBCT bite blocks or CBCT chin rests. Bag 'em baby! Our Bite Block Baggies are strong so they won't tear, and large enough to cover any bite block or chin rest you're using.

Economically priced too. Comes in a handy dispenser box.

500/box Item 70104

()

Comfee Ears

If your CBCT machine uses ear pieces to position the patient, you need Comfee Ears! These single use sleeves will fit the ear stabilization bars on most CBCT machines.

Give your patients the peace of mind they deserve.

150/box. Item 81500





۲



800.356.9729 | Fax: 631.242.1001 | www.flowdental.com





۲