


Hanes Mask – 3 Ply 100% Cotton with Antimicrobial Finish



COMPANY

America Knits is a textile manufacturer based in Swainsboro, Georgia. They are part of a consortium with Hanes to manufacture gowns and masks for healthcare providers as well as the federal government. Please note that Value First does not have a GPO agreement with this supplier; LeadingAge members are receiving the manufacturer cost with no fees being paid to Value First or any other entity outside of America Knits.

PRODUCT DESCRIPTION

Their mask has several properties that make it unique:

- Antimicrobial / antiviral properties impregnated into the fabric.
- The mask is reusable / washable. The antimicrobial properties are guaranteed to last 20 wash cycles. Mask have been washed 50 times and still retain their anti-microbial properties.
- The entire mask is made in the US. The raw material and production are all based in the US.
- The mask is FDA approved and is currently being provided to FEMA.
- Ask about adding your organization logo to your masks.

ORDERING INFORMATION

- The cost is \$1.75 per mask including shipping.
- A 50% deposit is required when the order is placed, and the balance is due 15 days after delivery of the product.
- Minimum order is 1 case, which is 500 masks.
- Complete the form on the next page to place an order.

It should be noted that these masks are not intended to be used when treating a resident that has a confirmed or suspected COVID-19 diagnosis. These masks are intended for general use by staff, residents and visitors that are not having direct interaction with individuals with confirmed cases of COVID-19. The masks are a great option based on the increased cost and lack of availability of disposable masks as well as the recent CMS guidance for everyone in the facility to wear a mask.

Mask Finished Measurements –

100% Cotton 4oz sq/yd Jersey Fabric with rubber inside cotton binding for stretch; mask is 3 ply over face
- pleats or fabric gathers on sides are acceptable if equipment restrictions exist

Finished Garment Measurements




All measurements are in inches				
POM #	Measurement Name	Tol (+)	Tol (-)	Spec (inches)
A	Length (at center)	0.500	0.500	5"
B	Length (at side)	0.500	0.500	3.5"
C	Width (at center)	0.500	0.500	7"
				7"
D	Ear Loop/Strap (total)	1.000	1.000	(7.25" if .0115 x 1/4" rubber is used)
				3.5"
DD	Ear Loop/Strap (half)	0.500	0.500	(3.625" if .0115 x 1/4" rubber is used)
				4.75"
DX	Ear Loop/Strap - Extended	0.500	0.500	(4.875" if .0115 x 1/4" rubber is used)
E	Binding Width/Height	0.125	0.125	0.375"
F	Distance between pleats	0.375	0.375	0.875"
G	Pleat Overlap	0.375	0.125	0.375"

*Measurements are Approximate

Garment has three pleats on right and left side (pleats are evenly spaced between top and bottom edge)
Note that Gathers are acceptable instead of pleats as long as fit is maintained
Right and Left side of garment should be symmetrical

Standard preferred rubber is HBI AK0649 = .0075" x 1/4"

**Approved substitute insert rubber is HBI AK0638 = .0115" x 1/4"



NOTE: Members should assess clinical acceptability of products prior to placing orders. Product availability changes regularly due to high demand, accessibility of raw materials and production capacity.